



Training on EQA and National TB Laboratory Network

Module 7: INTERPRETATION OF RECHECKING RESULTS AND FEEDBACK

Date

Uganda Supranational Reference Laboratory

Content Overview

- Feedback
- problem identification and solving

Feedback (1)

- Regular feedback
 - Indispensable
 - Encouragement and motivation
 - Improving performance, if applicable
 - No criticism unless serious negligence/bad will

Feedback (1) Cont....

- Visit by TB supervisor
 - Feedback; return of slides with serious errors
 - May identify possible causes of errors with checklist
 - May propose remedy, such as replacement of microscope

Feedback (2)

- Identify priority centres for problem solving
 - Several HFP and/or HFN
 - HFN and many LFN, or just many LFN
 - Many LFP
- On-site visit by an expert for identification of causes and problem solving

Problem Identification and Solving (1)

- In case of unsatisfactory performance of microscopy centres
 - Identify possible causes
 - Investigate
 - Take appropriate action
- Consider overall pattern of errors: (many) HFP, HFN, QE, combination of errors

Problem Identification and Solving (2)

- Many (or almost all) HFP and HFN
- Possible causes
 - Lack of training
 - Unusable microscope
 - Smears not examined
- Investigations
 - Examine 3+ with that microscope
 - Request staff to examine 2+, or 3+ and negative smear with good microscope

Problem Identification and Solving (3)

- Occasional HFP
- Possible causes
 - Administrative error
 - Failure of controllers
- Investigation
 - Compare lab. register with rechecking list for correct number and result

Problem Identification and Solving (4)

- Few HFP, with/without LFP
- Possible causes
 - Administrative error
 - Recognition of AFB?
 - No systematic restaining
 - Copying positive result 1st smear to 2nd and/or 3rd smear (workload?)

Problem Identification and Solving (4)

- Investigations

- Check admin. procedures: use of lab register, sputum request form, labeling sputum container
- Inconsistent smear results: regular isolated + or scanty?
- Restain and re-examine HFP to exclude that fading caused confusion and that it is not true FP

Problem Identification and Solving (5)

- Rare LFP
 - Ignore if at frequency comparable to other centres and controller

Problem Identification and Solving (6)

- Several LFP, with/without low grade HFP
- Possible causes
 - Bad microscope, artifacts taken as AFB?
 - Lack of experience
 - Contaminated carbol fuchsin stain
 - Poor counterchecks - FN by both controllers
- Investigations
 - Return FP to microscopy centre for showing AFB
 - Test carbol fuchsin on known negative smear

Problem Identification and Solving (7)

- Single HFN
- Possible causes
 - Administrative error
 - Smear not examined, 2nd and 3rd smear copied (workload?)
- Investigations
 - Check administrative procedures
 - Check workload
 - Smear not examined: difficult to prove

Problem Identification and Solving (8)

- Several HFN and/or LFN and/or low quantification
- Possible causes
 - Very thick smear
 - Poor light in microscope
 - Bad stain or poor staining (if smears were all restained)
 - Contaminated methylene blue or rinsing water (after restaining)
 - Superficial / no examination

Problem Identification and Solving (9)

- Several HFN and/or LFN and/or low quantification
- Investigations
 - Check thickness of smears
 - Check microscope: position condenser and diaphragm and remove filters
 - Check carbofuchsin stain: dark, red, shiny when poured on slide
 - Check if AFB are well stained in fresh not restained pos. smear
 - Check staining procedure: sufficient time, heating
 - Use methylene blue on known neg-atypical AFB?
 - Check workload (overload)
 - If all above fail to show cause; superficial reading may be the cause

Problem Identification and Solving (10)

- Serious QE
- Possible causes
 - Poor stain/staining
 - Problems with microscope, only few AFB in high positive smears
 - Lack of quantification skills
- Investigations
 - Check if AFB are well stained with carbol fuchsin
 - Observe staining procedure: sufficient time, heating
 - Ask staff to quantify a few positives

Summary

- Regular feedback to laboratories is very important for improving or maintaining motivation and performance
- A finding of poor performance has to be completed by a problem identification and solving visit. The problem and its cause(s) are rarely clear from rechecking as such.

Summary (Cont....)

- Feedback should include returning slides with discordant results to be reread by the peripheral laboratory personnel.
- All potential sources of error should be considered, including quality of stains and staining procedure, quality of microscopes, and administrative procedures that may contribute to recording errors.
- All problems contributing to errors must be resolved.

References

- WHO Laboratory Quality Management System Handbook
- WHO/GLI Tools.
- John, R. (1999). External Quality Assessment for AFB Smear Microscopy. *Public Health Practice Program Office Centers for Disease Control and Prevention, Rosemary Humes. Association of Public Health Laboratories, 17.*
- GLI Training package on EQA overview & Planning

Acknowledgments

