



Timely Accurate Diagonostics for a TB-Free Africa

### LABORATORY INFORMATION MANAGEMENT TRAINING

Module 4
Diagnostic connectivity solutions

SUPRANATIONAL REFERENCE LABORATORY, UGANDA

### **Module Objectives**

 Explain the concept of diagnostic connectivity and how diagnostic connectivity solutions contribute to epidemiological surveillance





### Diagnostic Connectivity Solutions

- Traditional diagnostics methods (e.g. smear microscopy and solid culture) typically require that a test result and associated patient details be manually written on paper, recorded in a register, and sent to the ordering clinician.
- \*\*Paper-based reporting is time consuming, and prone to transcription errors.





### Diagnostic Connectivity Solutions...

- \*Newer diagnostics (e.g GeneXpert®, Bactec MGIT) produce results data in digital format (also known as electronic data).
- Unlike written data on paper, electronic data can be rapidly and accurately sent to different recipients according to relevance reference Laboratory and utility and easily analysed.

Timely Accurate Diagonostics for a TB-Free Africa


#### **Diagnostic Connectivity** Solutions...

- Diagnostic Connectivity Solutions typically comprise of;
  - a connectable diagnostic device that produces electronic data
  - a software platform that receives and interprets data Supranational ®
    Reference Laboratory
  - a means to transmit data from




### Diagnostic Connectivity Solutions...

 Connectivity solutions facilitate the automatic transmission of electronic data to a variety of uses





### Diagnostic Connectivity Solutions





Timely Accurate Diagonostics for a TB-Free Africa

### Why Diagnostic Connectivity Solutions?

The We want to quickly and easily monitor and manage the performance of the machines, and report results







## Why Diagnostic Connectivity Solutions?...

- According to the WHO Framework of indicators and targets for laboratory strengthening under the End TB Strategy;
  - All sites that use WHO-recommended rapid diagnostics should be Supranational transmitting results electronically be required Discourse Discourse for a TR-Free Africa

### Why Diagnostic Connectivity Solutions?...

Remote monitoring via data connectivity solutions should be used to monitor key performance indicators at all sites that use WHOrecommended rapid diagnostics no later than 2020



- Remote monitoring and quality assurance
- Sending results automatically to clinicians
- Sending results to laboratory information management systems or electronic registers
- \*Inventory management





- Remote monitoring and quality assurance
  - Designated persons can use any internet-enabled computer to access the software platform, providing them with an overview of the facilities, devices and commodities in their network
    - For example: someone at a national reference laboratory can easily see how many tests are being performed and where, what are the results, and which sites are underperforming or experiencing abnormal results or errors



Reference Laboratory

ranational®

- Remote monitoring and quality assurance
- Sending results automatically to clinicians
  - Test results can automatically and instantly upon result availability be sent to a clinician's phone or email, SMS printer or other clinical results reporting mechanism, allowing for faster patient follow-up
    - A text message could also be sent to a patient, informing when their test results are ready and instructing them to visit the clinician to receive them





- Remote monitoring and quality assurance
- Sending results automatically to clinicians
- Sending results automatically to laboratory information management systems or electronic registers
  - Test results can be automatically integrated into laboratory information management systems or electronic registers, reducing staff time and the chance of transcription errors, and greatly facilitating M&E processes





- Remote monitoring and quality assurance
- Sending results automatically to clinicians
- Sending results automatically to laboratory information management systems or electronic registers
- Inventory management
  - Stocks can be entered at site-level and anticipated stock-out date or potential expiring cartridges can be forecast
  - Tracking of lot numbers can identify poor performance and abnormal error rates for quality assurance purposes





- · Remote monitoring and quality assurance
- Sending results automatically to clinicians
- Sending results automatically to laboratory information management systems or electronic registers
- Inventory management
- Surveillance
  - Real-time data and trends on disease or resistance patterns can be easily discerned using data from the network of devices






- Remote monitoring and quality assurance
- Sending results automatically to clinicians
- Sending results automatically to laboratory information management systems or electronic registers
- Inventory management
- Surveillance
- Data access
  - Subsets of data can be shared with partners or manufacturers



#### What is needed for Diagnostic Connectivity **Solutions** Software



- The selection of software depends on the needs and preferences of the Ministry of Health and may depend on:
  - functionalities of different systems
  - range of devices or languages supported
  - experience and capacity of the provider and implementing partners
  - data security measures and/or hosting arrangements in place
- Software may have complementary functionalities, justifying implementation and use of more than one system







Aspect

**Global Connectivity** 

Diagnostics

Platform (CDP)

# What is needed for Diagnostic Connectivity Solutions... Hardware

- Server
  - In-country, or 3<sup>rd</sup> party (cloud-based)
    - See slide on Hosting options, later in module
- Modem
  - Standalone desktop modems (also known as smart routers) are generally the best quality option but are more expensive
    - External antennas may provide stronger connection
  - USB internet dongles may be an option, but they carry the risk of being removed, misplaced and misused.





# What is needed for Diagnostic Connectivity Solutions Server

- There are two options for hosting of data collected:
  - Direct hosting using an in-country server
    - Requires adequate infrastructure, financial resources, and IT personnel to configure and maintain the server
  - 3<sup>rd</sup> party hosting (also known as: virtual or cloud-based servers)
    - Highly secure, provide redundancy against data loss, are quickly scalable and are a good option when a country does not have the IT capacity to properly configure and maintain an in-country server
    - Generally a cheaper option



### What is needed for **Diagnostic Connectivity** Solutions

- This Internet Connectivity
- The data plan (measured in Megabytes-MB or Gigabytes-GB) should be sufficient to allow the transmission of data from all connected devices to the server. The data plan should also be sufficient to Supranational® allow updates of anti-virus software



#### **Ensuring Impact from a Diagnostic Connectivity Solutions**

- Integration
  - Integration of diagnostics connectivity solutions with multiple types of diagnostics and with other databases (LIMS and electronic registers) allows Supranational® Supranational Supranation Su



#### **Ensuring Impact from a Diagnostic Connectivity** Solutions...

Sharing costs and responsibilities with other programmes = sustainability





#### Exercise

1. Explain the concept of diagnostic connectivity and how diagnostic connectivity solutions contribute to patient care and epidemiological surveillance





#### References

- GLI Quick Guide to TB Diagnostics
   Connectivity Solutions
- WHO/ERS. Digital health for the End TB Strategy: an agenda for action 2015
- WHO Guide on Electronic Recording and Reporting for Tuberculosis Care and



Control



#### **Acknowledgments**



















Timely Accurate Diagnostics for a TB-Free Africa



